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HYDROCARBONS, AROMATIC

1501

FORMULA: Table 1	MW: Table 1	CAS: Table 1	RTECS: Table 1
METHOD: 1501, Issue 3		EVALUATION: Full	
		Issue 1: 15 August 1990 Issue 3: 15 March 2003	
OSHA: Table 2 NIOSH: Table 2 ACGIH: Table 2		PROPERTIES: Table 1	
SYNONYMS: (Synonyms in Table 1)		Group A: benzene toluene ethylbenzene o-xylene m-xylene p-xylene Group B: cumene p-tert-butyltoluene α-methylstyrene β-methylstyrene styrene	
SAMPLING		MEASUREMENT	
SAMPLER:	SOLID SORBENT TUBE (coconut shell charcoal, 100 mg/50 mg)	TECHNIQUE:	GAS CHROMATOGRAPHY, FID
FLOW RATE:	Table 3	ANALYTE:	Hydrocarbons listed above
VOL-MIN:	Table 3	DESORPTION:	1 mL CS ₂ , stand 30 min with agitation
-MAX:	Table 3	INJECTION VOLUME:	1 μL (Group A: split 5:1; Group B: split 1:1)
SHIPMENT:	Routine	TEMPERATURE	
SAMPLE STABILITY:	30 days @ 5°C	-INJECTION:	250 °C
BLANKS:	10% of samples	-DETECTOR:	300 °C
		-COLUMN:	Group A: 40 °C (10 min) to 230°C (10 °C/min) Group B: 35°C (8 min) to 225°C (10°C/min)
		CARRIER GAS:	He @ 2.6 mL/min
ACCURACY		COLUMN:	Capillary, fused silica Group A: 30m x 0.32-mm ID; 1-μm film 100% PEG or equivalent Group B: 30m x 0.53-mm ID; 3-μm film crossbonded® 35% diphenyl 65% dimethyl polysiloxane or equivalent
RANGE STUDIED:	Table 3	CALIBRATION:	Solutions of analytes in CS ₂
BIAS:	Table 3	RANGE:	Table 4
OVERALL PRECISION (S _r):	Table 3	ESTIMATED LOD:	Table 4
ACCURACY:	Table 3	PRECISION (S _r):	Table 4
<p>APPLICABILITY: This method is for peak, ceiling, and TWA determinations of aromatic hydrocarbons. Interactions between analytes may reduce breakthrough volumes and affect desorption efficiencies. Naphthalene, originally validated in S292 [4], failed to meet acceptable desorption efficiency recovery and storage stability criteria at the levels evaluated in this study. However, the application of this method to naphthalene levels at or near the REL/PEL continues to meet acceptable recovery criteria. Styrene failed to meet acceptable recovery criteria at the two lowest levels evaluated in this study (highest level to meet the criteria was 181 μg/sample).</p>			
<p>INTERFERENCES: Under conditions of high humidity, the breakthrough volumes may be reduced. Other volatile organic compounds such as alcohols, ketones, ethers, and halogenated hydrocarbons are potential analytical interferences.</p>			
<p>OTHER METHODS: This method updates NMAM 1501 issued on August 15, 1994 [1] which was based upon P&CAM 127 (benzene, styrene, toluene, and xylene) [2]; S22 (p-tert-butyltoluene) [3]; S23 (cumene) [3]; S29 (ethylbenzene) [3]; S26 (α-methylstyrene) [3]; S30 (styrene); S311 (benzene) [4]; S343 (toluene) [4]; and S318 (xylenes) [4].</p>			

REAGENTS:

1. Carbon disulfide*, low benzene, chromatographic quality.
2. Analytes, reagent grade.
3. Helium, prepurified and filtered.
4. Hydrogen, prepurified and filtered.
5. Air, prepurified and filtered.

* See SPECIAL PRECAUTIONS

EQUIPMENT:

1. Sampler: glass tube, 7 cm long, 6-mm OD, 4-mm ID, flame-sealed ends, containing two sections of activated coconut shell charcoal (front = 100 mg, back = 50 mg) separated by a 2-mm urethane foam plug. A silylated glass wool plug precedes the front section and a 3-mm urethane foam plug follows the back section. Tubes are commercially available.
2. Personal sampling pump, 0.01 to 1.0 L/min (Table 3), with flexible connecting tubing.
3. Gas chromatograph, FID, integrator, and columns (page 1501-1).
4. Autosampler vials, glass, 1.8 mL, with PTFE-lined caps.
5. Pipets, 1-mL, and pipet bulb.
6. Syringes, 10- μ L, 25- μ L, and 250- μ L.
7. Volumetric flasks, 10-mL.

SPECIAL PRECAUTIONS: Carbon disulfide is toxic and extremely flammable (flash point = -30°C), benzene is a suspect carcinogen. Prepare standards and samples in a well ventilated hood.

SAMPLING:

1. Calibrate each personal sampling pump with a representative sampler in line.
2. Break the ends of the sampler immediately before sampling. Attach sampler to personal sampling pump with flexible tubing.
3. Sample at an accurately known flow rate between 0.01 and 0.2 L/min for a total sample size as shown in Table 3.
4. Cap the samplers with plastic (not rubber) caps and pack securely for shipment.

SAMPLE PREPARATION:

5. Place the front and back sorbent sections of the sampler tube in separate vials. Include the glass wool plug in the vial along with the front sorbent section.
6. Add 1.0 mL eluent to each vial. Attach crimp cap to each vial immediately.
7. Allow to stand at least 30 min with occasional agitation.

CALIBRATION AND QUALITY CONTROL:

8. Calibrate daily with at least six working standards from below the LOD to 10 times the LOQ. If necessary, additional standards may be added to extend the calibration curve.
 - a. Add known amounts of analytes to carbon disulfide solvent in 10-mL volumetric flasks and dilute to the mark. Prepare additional standards by serial dilution in 10-mL volumetric flasks.
 - b. Analyze together with samples and blanks (steps 11 through 12).
 - c. Prepare calibration graph (peak area of analyte vs. μ g analyte per sample).